



REF EU-24802-CVT (IPNXXXXXX)

8 Fr. 2 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

Arrowg+ard Blue® Two-Lumen CVC

Contents:

- 1: Arrowg+ard Blue® Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 30 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: ECG Extension Cable
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11


All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
 Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).
 Fluid path components are non-pyrogenic.
Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (14 Ga.)	0.88	4275
Proximal (14 Ga.)	0.90	2905

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4

LOT Lot Number
 YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Product of xxx
 Packaged in xxx

EU Authorized Representative and Importer: **Teleflex Medical IDA Business and Technology Park**, Dublin Road, Athlone Co. Westmeath, Ireland
CE 2797
 All components are CE 2797 unless otherwise noted.
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
Teleflex®
www.teleflex.com
 LBL059422 R02 (2024-10)



Arrowg+ard Blue® Two-Lumen CVC



REF EU-24802-CVT

8 Fr. 2 Lumen 30 cm catheter length .032 inch dia. spring-wire guide