



REF

CV-15122-UF

(IPNXXXXXX)

12<sup>Fr.</sup>

2<sup>Lumen</sup>

20<sup>cm catheter length</sup>

.035<sup>inch dia. spring-wire guide</sup>

# You-Bend™ Two-Lumen Hemodialysis Catheter for High Volume Infusions

- Contents:**
- 1: Two-Lumen Catheter†: 12 Fr. (4.06 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

‡MR Conditional

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

**Contraindications:** The large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.

- enbgcsdanelfi frde el hu it lt no
- plptro ru sk sl es sv tr

Lumen	Priming Volume* (mL)	Flow Rate (mL/min)	Venous Pressure Setting†† (mmHg)
Distal (12 Ga.)	1.8	200	100
Proximal (12 Ga.)	1.7	340	200
		200	100
		340	200

\* Priming volumes are approximate and are done without accessories.

†† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

EU Authorized Representative and Importer:

Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)  
LBL059702 R03 (2024-10)

Lot Number

YYYY-MM-DD

YYYY-MM-DD

(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number