



4^{Fr.} 2 Lumen

8 cm catheter length

.018 inch dia.
spring-wire
guide

Contents:

- 1: Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen


Priming
Volume*
(mL)Gravity
Flow
Rate†
(mL/hr)

Distal (22 Ga.)	0.28	548
Proximal (22 Ga.)	0.32	880

* Priming volumes are approximate and are done without accessories

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

en be bg hr cs da nl et fi fr de el hu is
it lv lt no pl pt ro ru sr sk sl es sv tr uk



LOT




EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060220 R01 (2023-11)



REF CS-15402-E

4 Fr. **2** Lumen

8 cm catheter length

.018 inch dia.
spring-wire
guide

ARROW

Pediatric Two-Lumen CVC