



REF CS-15402-E (IPNXXXXXX)

4 Fr. 2 Lumen 8 cm catheter length .018 inch dia. spring-wire guide

Pediatric Two-Lumen CVC

Contents:

- 1: Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm
 - 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
 - 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
 - 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
 - 1: Syringe: 5 mL Luer-Slip
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
 - 2: Dust Cap: Non-Vented
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
 - 1: SecondSite™ Adjustable Hub: Fastener
- All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.



Lumen


Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Distal (22 Ga.)	0.28	548
Proximal (22 Ga.)	0.32	880

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

en be bg hr cs da nl et fi fr de el hu is
it lv lt no pl pt ro ru sr sk sl es sv tr uk



All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Product of xxx
Packaged in xxx

teleflex®
www.teleflex.com

LBL060223 R02 (2024-08)



REF CS-15402-E

4 Fr. 2 Lumen

8 cm catheter length

.018 inch dia. spring-wire guide



Pediatric Two-Lumen CVC