



REF AW-04235

.035 inch dia. spring-wire guide

45 cm spring-wire guide

Marked Spring-Wire Guide with Arrow Advancer

Contents:

1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU

- en
- cs
- da
- nl
- et
- fi
- fr
- de
- el
- hu
- it
- no
- pl
- pt
- ro
- sk
- sl
- es
- sv



LOT



EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060574 R02 (2023-04)



REF AW-04235

.035 inch dia. spring-wire guide

45 cm spring-wire guide

Marked Spring-Wire Guide with Arrow Advancer

Contents:

1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU

- en
- cs
- da
- nl
- et
- fi
- fr
- de
- el
- hu
- it
- no
- pl
- pt
- ro
- sk
- sl
- es
- sv



LOT



EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060574 R02 (2023-04)