



REF AW-16402

.018 inch dia. spring-wire guide

68 cm spring-wire guide

Spring-Wire Guide with Arrow Advancer

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU

- en cs da nl et fi fr de el
- hu it no pl pt ro sk sl es
- sv



LOT



EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060630 R02 (2023-04)



REF AW-16402

.018 inch dia. spring-wire guide

68 cm spring-wire guide

Spring-Wire Guide with Arrow Advancer

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU

- en cs da nl et fi fr de el
- hu it no pl pt ro sk sl es
- sv



LOT



EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060630 R02 (2023-04)