

Product of xxx
Packaged in xxx



REF **SI-11142-S** (IPNXXXXXX)

9 Fr. distal lumen **2** Lumen **10** cm device length **.035** inch dia. spring-wire guide

MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

Contents:

- 1: Two-Lumen Central Venous Access Device: 4" (10 cm) Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 20.7 cm
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
- 1: Dust Cap: Vented
- 2: Dust Cap: Non-Vented
- 1: Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
 Fluid path components are non-pyrogenic.



Lumen

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.65	31440
Distal (9 Fr.) with 8 Fr. catheter	0.9	9040
Proximal (12 Ga.) no catheter	0.6	8970

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

en da fr de it pt es sv



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

EU Authorized Representative and Importer:
 Teleflex Medical IDA Business and Technology Park
 Dublin Road, Athlone Co. Westmeath, Ireland

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL060754 R03 (2024-03)



MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters



REF **SI-11142-S**

9 Fr. distal lumen **2** Lumen **10** cm device length **.035** inch dia. spring-wire guide