



REF PW-18130-P

.018 inch dia.
spring-wire
guide

130 cm
spring-wire
guide

Spring-Wire Guide .018" (0.46 mm) diameter x 51-3/16" (130 cm) Marked Wire

Contents:

1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex[®]
www.teleflex.com

LBL060986 R00 (2022-03)



REF PW-18130-P

.018 inch dia.
spring-wire
guide

130 cm
spring-wire
guide