



REF PW-18130

.018 inch dia.  
spring-wire  
guide

130 cm  
spring-wire  
guide

Spring-Wire Guide .018" (.46 mm) Diameter x 51-3/16" (130 cm) Wire. Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use.  
Failure to do so may result in severe patient injury or death.  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
www.teleflex.com

LBL060987 R00 (2022-03)



REF PW-18130

.018 inch dia.  
spring-wire  
guide

130 cm  
spring-wire  
guide



Spring-Wire Guide .018" (.46 mm) Diameter x 51-3/16" (130 cm) Wire. Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm