



REF CI-09800

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

**Contents:**

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.

All components are CE 2797 unless otherwise noted.

**Rx only**  
 Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
 Caution: This Product Contains Dry Natural Rubber.  
 Fluid path components are non-pyrogenic.

- en
- bg
- cs
- nl
- fr
- de
- el
- hu
- it
- pl
- pt
- ro
- ru
- sk
- sl
- es
- sv
- tr



EU Authorized Representative and Importer:  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
 LBL061033 R00 (2021-07)



REF CI-09800

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters