



REF CI-09800 (IPNXXXXXX)



8.5^{Fr.} 10^{cm} sheath length .035^{inch dia.} spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.

All components are CE 2797 unless otherwise noted.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Caution: This Product Contains Dry Natural Rubber.
Fluid path components are non-pyrogenic.


en bg cs nl fr de el hu it pl pt ro ru sk
sl es sv tr

EU Authorized Representative and Importer:
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.



LOT Lot Number
YYYY-MM-DD YYYY-MM-DD



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com
LBL061033 R01 (2025-06)



REF CI-09800

8.5^{Fr.} 10^{cm} sheath length .035^{inch dia.} spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/ Side Port for use with 7 - 7.5 Fr. Catheters

