



REF AW-05003

(IPNXXXXXX)

.032 inch dia. spring-wire guide

68 cm spring-wire guide

# Spring-Wire Guide

**Contents:**

- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark

[www.teleflex.com/IFU](http://www.teleflex.com/IFU)



**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)Lot Number

Product of xxx  
Packaged in xxx



EU Authorized Representative and Importer:



Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
 2797



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex**<sup>®</sup>  
[www.teleflex.com](http://www.teleflex.com)

LBL061106 R02 (2024-10)