



REF AW-05000

.032 inch dia. spring-wire guide

60 cm spring-wire guide

# Spring-Wire Guide

## Contents:


- 
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 

[www.teleflex.com/IFU](http://www.teleflex.com/IFU)



 EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
 2797

 Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL061107 R01 (2023-04)