



REF SI-09875-E

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
  - 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Pressure Transduction Probe
  - 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
  - 1: Obturator: 8 Fr.
  - 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
  - 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- All components are CE 2797 unless otherwise noted.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Caution: This Product Contains Dry Natural Rubber.  
Fluid path components are non-pyrogenic.

- en bg cs da nl fi fr de el hu it no pl pt
- ro ru sk sl es sv tr



EU Authorized Representative and Importer:  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL061116 R00 (2021-07)



REF SI-09875-E

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters