

Dimensions:
 Width: 8-3/4" (22,2 cm)
 Length: 11-1/8" (28,3 cm)

PB-10001-000A
PBZ-10001-000A

T-12703-100C
 Width: 8.7598 in = 222.5mm
 Length: 10.374 in = 263.5mm

ARROW



REF CS-17752 (IPNXXXXXX)

7 Fr. **2** Lumen **60** cm catheter length **.035** inch dia. spring-wire guide

Two-Lumen Antecubital CVC

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 60 cm
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

CAS # 7440-48-4

en be bg hr cs da nl et fi fr de el hu is it lv lt no pl pt ro ru sr
 sk sl es sv tr uk

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (14 Ga.)	1.13	2216
Proximal (18 Ga.)	0.80	396

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

EC/REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
CE 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Product of xxx
 Packaged in xxx

Teleflex
www.teleflex.com

LBL061125 R02 (2025-01)



Two-Lumen Antecubital CVC

ARROW

REF CS-17752

7 Fr. **2** Lumen **60** cm catheter length **.035** inch dia. spring-wire guide