



REF AW-06832-N

(IPNXXXXXX)

**.032** inch dia.  
spring-wire  
guide

**68** cm  
spring-wire  
guide

# Spring-Wire Guide

Contents: 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG mark

[www.teleflex.com/IFU](http://www.teleflex.com/IFU)

en cs da nl et fi fr de el hu it  
no pl pt ro sk sl es sv



EU Authorized Representative and Importer:  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

LOT Lot Number

YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)

LBL061132 R02 (2024-10)