



REF AW-06032-N

(IPNXXXXXX)

.032 inch dia.
spring-wire
guide

60 cm
spring-wire
guide

Spring-Wire Guide

Contents: 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark

www.teleflex.com/IFU

en cs da nl et fi fr de el hu it
no pl pt ro sk sl es sv



LOT Lot Number

YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx
Packaged in xxx



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL061133 R02 (2024-10)