

Dimensions:  
 Width: 6-9/64" (15.6 cm)  
 Length: 8-21/32" (22 cm)

**PB-50014-100B**  
**PBZ-50014-103A**

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



**REF** CV-14702  
 (IPNXXXXXX)

**ARROW**

# Two-Lumen CVC

**Contents:**


- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.** Fluid path components are non-pyrogenic.



CAS # 7440-48-4

EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Product of xxx  
 Packaged in xxx

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.57	1917
Proximal (16 Ga.)	0.62	2000

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)Lot Number

LBL061184 R04 (2024-10)



**7** Fr. **2** Lumen **30** cm catheter length **.032** inch dia. spring-wire guide

**REF** CV-14702

**Two-Lumen CVC**

**ARROW**