



REF HO-14702 (IPNXXXXXX)

7 Fr. 2 Lumen

30 cm catheter length












23 cm insertable catheter length with Cath-Gard®


HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield origin:xxxxxxxxxxxxxxxxxxxxx
2: Dust Cap: Non-Vented origin:xxxxxxxxxxxxxxxxxxxxx
All components are CE 2797 unless otherwise noted.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.



CAS # 7440-48-4

en

be

bg

hr

cs

da

nl

et

fi

fr

de

el

hu

is

it

lv

lt

no

pl

pt

ro

ru

sr

sk


sl

es

sv

tr

uk


 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.57	1917
Proximal (16 Ga.)	0.62	2000

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Teleflex®
www.teleflex.com

LBL061678 R02 (2024-05)



7 Fr. 2 Lumen

30 cm catheter length

23 cm insertable catheter length with Cath-Gard®

REF HO-14702

HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System

