



REF HO-14703

7 Fr. 3 Lumen

30 cm catheter length

23 cm insertable catheter length with Cath-Gard®

**HANDS-OFF® Three-Lumen CVC for use only with Arrow One-Piece PSI System**

**Contents:**

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield
- 3: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**  
**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.



CAS # 7440-48-4



| Lumen             | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (16 Ga.)   | 0.54                 | 2025                       |
| Medial (18 Ga.)   | 0.45                 | 896                        |
| Proximal (18 Ga.) | 0.48                 | 955                        |

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:



Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL061693 R01 (2023-08)

LOT



7 Fr. 3 Lumen

30 cm catheter length

23 cm insertable catheter length with Cath-Gard®

REF HO-14703

**HANDS-OFF® Three-Lumen CVC for use only with Arrow One-Piece PSI System**