

ARROW®



REF EU-02041-ML (IPNXXXXXX)

4 Fr. | 1 Lumen | 20 cm catheter length | .018 inch dia. spring-wire guide

# Arrow® Midline Kit


**Contents:**

- 1: Midline Catheter, Trimmable: 4 Fr. x 20 cm Radiopaque Polyurethane with T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire
- 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Lock
- 1: SharpsAway® II Locking Disposal Cup
- 1: Catheter Trimmer
- 1: Safety Scalpel: #11
- 1: Paper Tape Measure
- 1: Tourniquet
- 1: Dressing: GRIP-LOK®<sup>1</sup> Catheter Fixation Device
- 1: Valve: MicroCLAVE®<sup>2</sup> Neutral Displacement Connector

<sup>1</sup>A registered trademark of TIDI.  
<sup>2</sup>A registered trademark of ICU Medical, Inc.

- nl** Arrow midline katheterkit
- fr** Kit Arrow Midline
- de** Arrow Midline-Kit
- it** Kit catetere Midline Arrow
- pt** Kit de Linha Mediana da Arrow
- es** Kit de catéter de línea media de Arrow
- sv** Arrow Midline-sats

Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.30	1690	7110

\* Priming volumes are approximate and are done without the injection cap.  
† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.  
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
Above values apply to the untrimmed catheter. Trimming catheter will alter priming volumes and flow rates.



LOT LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
www.teleflex.com  
LBL062089 (2021-06)



REF EU-02041-ML

4 Fr. | 1 Lumen | 20 cm catheter length | .018 inch dia. spring-wire guide

# Arrow® Midline Kit

ARROW®