

**DO NOT BEND****ARROW**REF **SAC-01618-PBX** (IPNXXXXXX)**18** Ga.**16** cm catheter length**Arterial Catheterization Set****Contents:**

- 1: Indwelling Catheter: 18 Ga. x 16 cm Radiopaque PEBA  
 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer  
 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

Catheter Priming Volume: 0.24 mL Priming Volume is approximate.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- cs** Arteriální katetrizační set
- da** Arterielt kateteriseringssæt
- nl** Arteriële katheterisatieset
- et** Arterikanüüli komplekt
- fi** Valtimokatetrisarja
- fr** Jeu de cathétérisme artériel
- de** Besteck zur Katheterisierung von Arterien
- el** Σετ αρτηριακού καθετηριασμού
- hu** Artériás katéter készlet
- it** Set per cateterismo arterioso
- lt** Arterijų kateterizavimo rinkinys
- no** Arteriell kateteriseringssett
- pl** Zestaw do cewnikowania tętnic
- pt** Conjunto de cateterização arterial
- ro** Set de cateterizare arterială
- ru** Набор для артериальной катетеризации
- sk** Arteriálna katetrizačná sada
- sl** Pribor za arterijsko kateterizacijo
- es** Conjunto de cateterización arterial
- sv** Artärkateteriseringssset
- tr** Arteriyel Kateterizasyon Seti



(01)XXXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (10)LotNumber



LotNumber



YYYY-MM-DD

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex**

www.teleflex.com

LBL062149 (2021-06)



Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

REF **SAC-01618-PBX****ARROW****Arterial Catheterization Set****18** Ga.**16** cm catheter length**DO NOT BEND**