



REF FA-04020

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 9-5/8" (24 cm)
- 1: Catheter: 20 Ga. x 4-1/4" (10.8 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Wing Clip: Suture

Rx only



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA



teleflex®
 www.teleflex.com

ARROW

LBL062856 R00 (2022-01)



REF FA-04020

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 9-5/8" (24 cm)
- 1: Catheter: 20 Ga. x 4-1/4" (10.8 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Wing Clip: Suture

Rx only



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA



teleflex®
 www.teleflex.com

ARROW

LBL062856 R00 (2022-01)



REF FA-04020

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 9-5/8" (24 cm)
- 1: Catheter: 20 Ga. x 4-1/4" (10.8 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Wing Clip: Suture

Rx only



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA



teleflex®
 www.teleflex.com

ARROW

LBL062856 R00 (2022-01)



REF FA-04020

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .018" (0.46 mm) dia. x 9-5/8" (24 cm)
- 1: Catheter: 20 Ga. x 4-1/4" (10.8 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Wing Clip: Suture

Rx only



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA



teleflex®
 www.teleflex.com

ARROW

LBL062856 R00 (2022-01)