

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)

ARROW®



REF GH-04120

(IPNXXXXXX)

20^{Ga.} 6.35^{cm} .021^{inch dia.}
catheter length spring-wire guide

Radial Artery Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 6.35 cm Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



CAS # 7440-48-4

EU Authorized
Representative and
Importer:Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland**All components are CE 2797 unless otherwise noted.**Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USAProduct of xxx
Packaged in xxx**Teleflex®**
www.teleflex.com

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL062860 R01 (2025-05)