



REF GH-04122

22<sup>Ga.</sup> 3.49<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# Radial Artery Catheterization Set

## Contents:

- 1: Catheter: 22 Ga. x 3.49 cm Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW

All components are CE 2797 unless otherwise noted.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) Not made with natural rubber latex. Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL062862 R00 (2022-01)

LOT

