



REF SI-09600

6 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 4 - 6 Fr. Catheters

**Contents:**

- 1: Arrow-Flex® Sheath: 6 Fr. (1.9 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 30.5 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)


All components are CE 2797 unless otherwise noted.

**Rx only**

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Caution: This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic.



LOT



EU Authorized Representative and Importer:

Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex**  
www.teleflex.com

LBL063945 R00 (2022-02)



REF SI-09600

6 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 4 - 6 Fr. Catheters