



REF SI-09700

7 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 Fr. Catheters

**Contents:**

- 1: Arrow-Flex® Sheath: 7 Fr. (2.3 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 30.5 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.

**Rx only**

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Caution: This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL063950 R00 (2022-02)



REF SI-09700

7 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 Fr. Catheters