



REF

SI-09806

(IPNXXXXXX)



8.5^{Fr.}

10^{cm}
sheath
length

.025<sup>inch dia.
spring-wire
guide</sup>

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Caution: This Product Contains Dry Natural Rubber.
Fluid path components are non-pyrogenic.

enbgcsnlfifrdedeluhuitltonoplpt

roaruskslesvstr

EU Authorized Representative and Importer:
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

EC REP

CE 2797

All components are CE 2797 unless otherwise noted.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx.
Packaged in xxx.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com
LBL063954 R01 (2024-10)



REF

SI-09806

8.5^{Fr.}

10^{cm}
sheath
length

.025<sup>inch dia.
spring-wire
guide</sup>

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/ Side Port for use with 7 - 7.5 Fr. Catheters

