

REF

RC-05801

(IPNXXXXXX)



Trauma Kit

8.5^{Fr.}

8.89^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

Contents:

- 1: Sheath: 8.5 Fr. x 3-1/2" (8.89 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Syringe: 30 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: OpSite®1 5-1/2" x 4" (14 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of T.J. Smith & Nephew, Limited.
All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.





EU Authorized
Representative and
Importer:


EC REP

CE

2797

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.








Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx.
Packaged in xxx






Teleflex®
www.teleflex.com

LBL064064 R02 (2024-10)






STERILE




LOT


Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

REF

RC-05801

8.5^{Fr.}

8.89^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

Trauma Kit

