

REF RC-05801

Trauma Kit

8.5^{Fr.} 8.89^{cm} sheath length .035^{inch dia.} spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 3-1/2" (8.89 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Syringe: 30 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 1: Dressing: OpSite®1 5-1/2" x 4" (14 cm x 10 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of T.J. Smith & Nephew, Limited.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



LOT



EU Authorized Representative and Importer:

EC REP
 CE 2797
 Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL064064 R00 (2022-02)

REF RC-05801

Trauma Kit

8.5^{Fr.} 8.89^{cm} sheath length .035^{inch dia.} spring-wire guide