

REF ST-09875

ARROW[®]

TwistLock™ Cath-Gard® Catheter Contamination Shield for use with 7.5 and 8 Fr. Catheters



Contents:

1: Cath-Gard Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

Not made with natural rubber latex. Fluid path components are non-pyrogenic.

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex[®]
www.teleflex.com

 EU Authorized
Representative and
Importer:

LBL064183 R00 (2022-03)

