



REF FA-04018 (IPNXXXXXX)

Rx only

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 9-21/32" (24.5 cm)
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 4-1/4" (10.8 cm) Radiopaque over 20 Ga. TW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture_origin:xxxxxxxxxxxxxxxxxxxxxxxx

CAS # 7440-48-4 Set Product of xxx
 Packaged in xxx
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
 www.teleflex.com

LOT Lot Number
 YYY-YY-MM-DD
 YYY-YY-MM-DD

CE 2797



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL067280 R02 (2024-05)



REF FA-04018 (IPNXXXXXX)

Rx only

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 9-21/32" (24.5 cm)
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 4-1/4" (10.8 cm) Radiopaque over 20 Ga. TW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture_origin:xxxxxxxxxxxxxxxxxxxxxxxx

CAS # 7440-48-4 Set Product of xxx
 Packaged in xxx
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
 www.teleflex.com

LOT Lot Number
 YYY-YY-MM-DD
 YYY-YY-MM-DD

CE 2797



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL067280 R02 (2024-05)



REF FA-04018 (IPNXXXXXX)

Rx only

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 9-21/32" (24.5 cm)
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 4-1/4" (10.8 cm) Radiopaque over 20 Ga. TW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture_origin:xxxxxxxxxxxxxxxxxxxxxxxx

CAS # 7440-48-4 Set Product of xxx
 Packaged in xxx
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
 www.teleflex.com

LOT Lot Number
 YYY-YY-MM-DD
 YYY-YY-MM-DD

CE 2797



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL067280 R02 (2024-05)



REF FA-04018 (IPNXXXXXX)

Rx only

Arterial Catheterization Set

Contents:

- 1: Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 9-21/32" (24.5 cm)
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 4-1/4" (10.8 cm) Radiopaque over 20 Ga. TW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture_origin:xxxxxxxxxxxxxxxxxxxxxxxx

CAS # 7440-48-4 Set Product of xxx
 Packaged in xxx
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
 www.teleflex.com

LOT Lot Number
 YYY-YY-MM-DD
 YYY-YY-MM-DD

CE 2797



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL067280 R02 (2024-05)