

Dimensions:  
Width: 8-3/4" (22,2 cm)  
Length: 11-1/8" (28,3 cm)

PB-10001-000A  
PBZ-10001-000A

T-12703-100C  
Width: 8.7598 in = 222.5mm  
Length: 10.374 in = 263.5mm

ARROW



REF PW-18080

(IPNXXXXXX)

.018 inch dia.  
spring-wire  
guide

80 cm  
spring-wire  
guide

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Fluid path components are non-pyrogenic.



EU Authorized  
Representative and  
Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Product of xxx  
Packaged in xxx

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL068225 R01 (2024-07)



REF PW-18080

.018 inch dia.  
spring-wire  
guide

80 cm  
spring-wire  
guide

ARROW

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm