



REF PW-18080

.018 inch dia.
spring-wire
guide

80 cm
spring-wire
guide

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
 Not made with natural rubber latex.
 Fluid path components are non-pyrogenic.



LOT



EU Authorized Representative and Importer:

Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 2797

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL068225 R00 (2023-02)



REF PW-18080

.018 inch dia.
spring-wire
guide

80 cm
spring-wire
guide

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm