



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)



REF GH-04120-E

LOT



### Radial Artery Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) Radiopaque Polyurethane with Integral Suture Wing
  - 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 7-1/2" (19 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Rx only



**Teleflex**<sup>®</sup>  
 www.teleflex.com  
 LBL068459 R00 (2022-09)

