

Product of xxx  
Packaged in xxx



REF **EU-25063-HPMSB** (IPNXXXXXX)

**6** Fr. **3** Lumen **50** cm catheter length **.018** inch dia. spring-wire guide

# Pressure Injectable Three-Lumen PICC

**Contents:**

- 1: Three Lumen Catheter: 6 Fr. (2.14 mm OD) x 50 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire
- 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
- 1: Injection Needle: Eclipse™<sup>1</sup> 22 Ga. x 1-5/8" (4 cm)
- 1: Injection Needle: Eclipse™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
- 2: Syringe: 10 mL Luer-Lock
- 1: Syringe: 3 mL Luer-Lock
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Catheter Trimmer
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 60" x 76" (152 cm x 193 cm)
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel
- 1: Safety Scalpel: #11
- 1: Patient ID Card
- 1: Sterile Procedure Sign
- 1: Patient Information Booklet
- 1: Chart Sticker
- 1: PICC Line Insertion Checklist
- 1: Paper Tape Measure
- 1: Tourniquet
- 5: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 10: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®<sup>2</sup> I.V. Advanced 3-1/4" x 4-1/2" (8.5 cm x 11.5 cm)
- 1: Dressing: GRIP-LOK®<sup>3</sup> Catheter Fixation Device
- 1: Surgical Apparel: Cap
- 1: Surgical Apparel: Impervious Gown
- 1: Dressing: Steri-Strip®<sup>2</sup> 1/2" x 4" (1.2 cm x 10 cm)
- 1: Surgical Apparel: Mask
- 1: Surgical Apparel: Mask with Shield
- 1: Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)
- 3: Valve: MicroCLAVE™<sup>4</sup> Neutral Displacement Connector
- 1: HemoHopper® Fluid Receptacle

<sup>2</sup>A registered trademark of 3M Company.  
<sup>3</sup>A registered trademark of TIDI.  
<sup>4</sup>A trademark of ICU Medical, Inc.  
**All components are CE 2797 unless otherwise noted.**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

**Contraindications:** The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (17 Ga.)	0.62	1195	6410	6
Medial (19 Ga.)	0.42	297	2720	not rated
Proximal (19 Ga.)	0.45	297	2830	not rated

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.  
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
\*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



CAS # 7440-48-4

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:  
**Teleflex Medical**  
IDA Business and Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
LBL069885 R02 (2024-02)



# Pressure Injectable Three-Lumen PICC

REF **EU-25063-HPMSB**

**6** Fr. **3** Lumen **50** cm catheter length **.018** inch dia. spring-wire guide

