

Dimensions:  
Width: 6-9/64" (15,6 cm)  
Length: 8-21/32" (22 cm)

PB-50014-100B  
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

Product of xxx  
Packaged in xxx

ARROW®



REF CS-15403 (IPNXXXXXX)

Pediatric Three-Lumen CVC

Contents:



- 1: Three-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen

Priming  
Volume\*  
(mL)

Gravity  
Flow  
Rate†  
(mL/hr)

Distal (20 Ga.)	0.28	601
Medial (23 Ga.)	0.30	106
Proximal (23 Ga.)	0.34	119

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4



EU Authorized  
Representative and  
Importer:



2797

Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL074338 R00 (2024-01)

LOT LotNumber



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber



4 Fr. 3 Lumen

8 cm  
catheter  
length

.018 inch dia.  
spring-wire  
guide

REF CS-15403

Pediatric Three-Lumen CVC

ARROW®