



AGB⁺

REF

EU-42802-HPSKM
(IPNXXXXXX)

8 Fr. 2 Lumen

16 cm catheter length

.032 inch dia. spring-wire guide

ARROW[®]

Pressure Injectable Arrowg+ard Blue Plus[®] Two-Lumen CVC

Contents:

- 1: Arrowg+ard Blue Plus[®] 1 Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 16 cm, Pressure Injectable

1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow[®] GlideWheel[™] Wire Advancer with ECG Mark

1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

1: Injection Needle: Eclipse[™] 25 Ga. x 1" (2.54 cm)

1: Pressure Transduction Probe

1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip

1: Syringe: 3 mL Luer-Lock

1: Syringe: 10 mL Luer-Slip

1: Arrow[®] Raulerson Spring-Wire Introduction Syringe: 5 mL

1: Tissue Dilator: 9 Fr. x 10.2 cm

2: Dust Cap: Non-Vented

1: SecondSite[™] Adjustable Hub: Fastener

1: SecondSite[™] Adjustable Hub: Catheter Clamp

1: SharpsAway[®] II Locking Disposal Cup

1: Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation

1: Towel

1: Safety Scalpel: #11

3: Prep Sponge Swab

10: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Dressing: Tegaderm[®] 3 I.V. Advanced 3-1/4" x 4-1/2" (8.5 cm x 11.5 cm)

1: Dressing: GRIP-LOK[®] Catheter Fixation Device

1: Surgical Apparel: Cap
- 1: Surgical Apparel: Impervious Gown

1: Surgical Apparel: Mask with Shield

1: Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)

1: Tray: Prep

1: Licensed under US Patent No. 7,329,412.

2: A trademark of Becton, Dickinson and Company.

3: A registered trademark of 3M Company.

4: A registered trademark of TID.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.


Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).
Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (14 Ga.)	0.68	6350	10
Proximal (14 Ga.)	0.70	4980	10

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.
** Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.



LOT Lot Number
YYYY-MM-DD
YYYY-MM-DD



01XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

CAS # 7440-48-4



EU Authorized Representative and Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex[®]
www.teleflex.com

LBL075232 R01 (2026-01)



REF

EU-42802-HPSKM

LOT

Lot Number

(IPNXXXXXX)



YYYY-MM-DD



YYYY-MM-DD



REF

EU-42802-HPSKM

LOT

Lot Number

(IPNXXXXXX)



YYYY-MM-DD



YYYY-MM-DD