

Product of xxx
Packaged in xxx

ARROW



REF AW-04032

(IPNXXXXXX)

DuoFlex® Spring-Wire Guide Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide: .032" (0.81 mm) dia. x 12" (30 cm) (Straight Soft Tip on Both Ends)

Rx only www.teleflex.com/IFU



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL076480 R00 (2026-05)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx

ARROW



REF AW-04032

(IPNXXXXXX)

DuoFlex® Spring-Wire Guide Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide: .032" (0.81 mm) dia. x 12" (30 cm) (Straight Soft Tip on Both Ends)

Rx only www.teleflex.com/IFU



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL076480 R00 (2026-05)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx

ARROW



REF AW-04032


(IPNXXXXXX)

DuoFlex® Spring-Wire Guide Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide: .032" (0.81 mm) dia. x 12" (30 cm) (Straight Soft Tip on Both Ends)

Rx only www.teleflex.com/IFU



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL076480 R00 (2026-05)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx

ARROW



REF AW-04032

(IPNXXXXXX)

DuoFlex® Spring-Wire Guide Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide: .032" (0.81 mm) dia. x 12" (30 cm) (Straight Soft Tip on Both Ends)

Rx only www.teleflex.com/IFU



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL076480 R00 (2026-05)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD