

Product of xxx
Packaged in xxx

ARROW



REF AW-04432

(IPNXXXXXX)

.032 inch dia.
spring-wire
guide


60 cm
spring-wire
guide

Marked Spring-Wire Guide with Arrow Advancer Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J"
Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL076490 R00 (2026-05)

Product of xxx
Packaged in xxx

ARROW



REF AW-04432

(IPNXXXXXX)

.032 inch dia.
spring-wire
guide


60 cm
spring-wire
guide

Marked Spring-Wire Guide with Arrow Advancer Maximum OD per ISO 11070 is 0.82 mm.

Contents:

1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J"
Tip on Other) with Arrow Advancer

Rx only www.teleflex.com/IFU



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL076490 R00 (2026-05)