



REF HO-14703 (IPNXXXXXX)

7 Fr. 3 Lumen 30 cm catheter length 23 cm insertable catheter length with Cath-Gard®

HANDS-OFF® Three-Lumen CVC for use only with Arrow One-Piece PSI System

Contents:

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield
 - 3: Dust Cap: Non-Vented
- All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
 Not made with natural rubber latex.
 Fluid path components are non-pyrogenic.



CAS # 7440-48-4



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.54	2025
Medial (18 Ga.)	0.45	896
Proximal (18 Ga.)	0.48	955

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Product of xxx
 Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Teleflex®
 www.teleflex.com

LBL076627 R00 (2025-09)



7 Fr. 3 Lumen 30 cm catheter length 23 cm insertable catheter length with Cath-Gard®

REF HO-14703

HANDS-OFF® Three-Lumen CVC for use only with Arrow One-Piece PSI System

