



REF IC-09880

(IPNXXXXXX)



8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

PSI Set for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Dust Cap: Non-Vented
- 1: Safety Scalpel: #11
- 1: Stopcock: 3-way High-Flow
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
Caution: This Product Contains Dry Natural Rubber.
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL076628 R00 (2025-01)



REF IC-09880

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

PSI Set for use with 7 - 7.5 Fr. Catheters

