



REF

SI-09600

(IPNXXXXXX)



6^{Fr.}

10^{cm}
sheath
length

.035<sup>inch dia.
spring-wire
guide</sup>

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 4 - 6 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 6 Fr. (1.9 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 30.5 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Caution: This Product Contains Dry Natural Rubber.
Fluid path components are non-pyrogenic.

- enbgcsnlfifrdedeluhuitltonoplpt
- roaruskslesvstr

EU Authorized Representative and Importer:
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.



Lot Number

YYYY-MM-DD

YYYY-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx.
Packaged in xxx.

Teleflex®
www.teleflex.com
LBL076685 R01 (2025-07)



REF

SI-09600

6^{Fr.}

10^{cm}
sheath
length

.035<sup>inch dia.
spring-wire
guide</sup>

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/ Side Port for use with 4 - 6 Fr. Catheters

