



REF SI-09700

(IPNXXXXXX)

10 cm sheath length

.035 inch dia. spring-wire guide



Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 Fr. Catheters

- Arrow-Flex® Sheath: 7 Fr. (2.3 mm tip ID) x 10 cm and Tissue Dilator
- Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- Obturator: 8 Fr.
- Cath-Gard® Catheter Contamination Shield: 30.5 cm with Sterile Tape Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Caution: This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic





EU Authorized Representative and Importer:



Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA


Product of xxx Packaged in xxx

Teleflex®

LBL076687 R01 (2025-07)







STERILE EO



(17)YYMMDD (11)YYMMDD (10)Lot Number





(A) CAS # 7440-48-4

