



(IPNXXXXXX)





.035 inch dia. spring-wire guide

## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

## Contents:

- Arrow-Flex® Sheath: 8 Fr. (2.5 mm tip ID) x 10 cm and Tissue Dilator Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

- Obturator: 8 Fr.
  Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
  Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted

## Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe

patient injury or death. www.teleflex.com/IFU
Caution: This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic.























EU Authorized Representative and Importer:




Teleflex Medical IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Product of xxx Packaged in xxx







(01)XXXXXXXXXXXXXXX (11)YYMMDD (10)Lot Number

**Meleflex**®

LBL076689 R00 (2025-01)



035 inch dia. spring-wire guide

