

Dimensions:  
Width: 10-5/8"  
Length: 12-1/2"

PB-10003-000

T04400101A3 Dimensions: 10 1/2" x 10.5"

K08800109A3 Dimensions: 10 1/2" x 11"



REF SI-09880 (IPNXXXXXX)



8.5<sup>Fr.</sup> 10<sup>cm</sup> sheath length .035<sup>inch dia.</sup> spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

All components are CE 2797 unless otherwise noted.


Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

en bg cs nl fi fr de el hu it lt no pl pt  
ro ru sk sl es sv tr

EU Authorized Representative and Importer:  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



LOT Lot Number

YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Product of xxx  
Packaged in xxx

Teleflex®  
www.teleflex.com  
LBL076693 R00 (2025-01)



REF SI-09880

8.5<sup>Fr.</sup> 10<sup>cm</sup> sheath length .035<sup>inch dia.</sup> spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

