

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)

ARROW®

REF RC-09850

(IPNXXXXXX)

8.5^{Fr.}

6.4 cm catheter length

.025 inch dia.
spring-wire
guide

RIC® Rapid Infusion Catheter Exchange Set

Contents:

- 1: Sheath: 8.5 Fr. x 2-1/2" (6.4 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator
1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com



Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL077039 R00 (2025-09)