



REF **CS-15402-E** (IPNXXXXXX)

4 Fr. **2** Lumen **8** cm catheter length **.018** inch dia. spring-wire guide

Pediatric Two-Lumen CVC

Contents:

- 1: Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm
 - 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
 - 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
 - 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
 - 2: Dust Cap: Non-Vented
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
 - 1: SecondSite™ Adjustable Hub: Fastener
- All components are CE 2797 unless otherwise noted.**

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen

Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Distal (22 Ga.)	0.28	548
Proximal (22 Ga.)	0.32	880

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4

EU Authorized Representative and Importer:
EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone Co. Westmeath, Ireland
CE 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Product of xxx
 Packaged in xxx

Teleflex®
www.teleflex.com
 LBL077126 R00 (2025-07)



REF **CS-15402-E**

4 Fr. **2** Lumen

8 cm catheter length

.018 inch dia. spring-wire guide



Pediatric Two-Lumen CVC