



REF CS-22703-E

(IPNXXXXXX)



Arrowg+ard Blue® Three-Lumen CVC

Contents:

- 1: Arrowg+ard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
 - 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
 - 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
 - 1: Syringe: 5 mL Luer-Slip
 - 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
 - 3: Dust Cap: Non-Vented
 - 1: SecondSite™ Adjustable Hub: Fastener
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
- All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (16 Ga.) | 0.42 | 2656 |
| Medial (18 Ga.) | 0.38 | 1353 |
| Proximal (18 Ga.) | 0.41 | 1477 |

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer:

EC REP
 Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Product of xxx
 Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

LBL077135 R00 (2025-08)



7 Fr. 3 Lumen

16 cm catheter length

.032 inch dia. spring-wire guide

REF CS-22703-E

Arrowg+ard Blue® Three-Lumen CVC

