

Dimensions:
Width: 8-3/4" (22,2 cm)
Length: 11-1/8" (28,3 cm)

PB-10001-000A
PBZ-10001-000A

T-12703-100C
Width: 8.7598 in = 222.5mm
Length: 10.374 in = 263.5mm

ARROW



REF PW-18080

(IPNXXXXXX)

.018 inch dia.
spring-wire
guide


80 cm
spring-wire
guide

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Fluid path components are non-pyrogenic.



EU Authorized
Representative and
Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx
Packaged in xxx

Teleflex
www.teleflex.com

LBL077151 R00 (2025-09)



REF PW-18080

.018 inch dia.
spring-wire
guide

80 cm
spring-wire
guide

ARROW

Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm